

WEST PALM BEACH FIRE DEPARTMENT



BACTERIAL/VIRAL FILTERS



PURPOSE:

- Filter particulates. **SINGLE USE ONLY.**
- **The 19mm filters are intended for the CPAP.**
- **The 15mm/22mm filters are intended for BVM/ I-GEL/ETT.**

EQUIPMENT:

- One (1) disposable 19mm Bacterial/Viral Filter.
- One (1) disposable 15mm/22mm Bacterial/Viral Filter

DIRECTIONS OF USE:

- Place filter closest to the patient.
- Ensure connectors are secure and free from leaks.
- If utilizing a mask ensure a proper seal to reduce the transmission of microbes.



**Filters up to 99.99% of expiratory
bacterial/viral content**

WPBFD
EST. 1894